

**WASH ALL GEL
PRODUCTS BY
HAND WITH MILD
SOAP AND WATER**



AIR DRY

ALL SIZES ARE APPROXIMATE				
Product Number	Product Description	Circumferential Measurements		For a snugger fit, go down a size For a looser fit go up a size
		Relaxed (in)	Stretched (in)	
300P	Gel Toe/Finger Cap S/M	1 7/8	3	
301P	Gel Toe/Finger Cap L/XL	2 1/2	3 3/4	
302P	Gel Toe/Finger Cap XXL	3	4 1/2	
306P	Gel Tubing Narrow 3/4" x 6"	1 3/4	2 5/8	
307P	Gel Tubing Wide 1" x 6"	2 1/2	4 1/4	
	Gel Tubing Extra Wide	6"	10"	
K1025	Half Gel Tubing Small	5.5"	16"	
K1026	Half Gel Tubing Medium	7.5"	18"	
K1027	Half Gel Tubing Large	9.5"	20"	
1001P	Gel Ankle Sleeve 5" x 10"	6	10	
1008P	Ankle Gel Disk Sleeve	6 3/4	10	measured at Ankle
1004P	Gel Heel Sleeve S/M	6 3/4	10	measured at Ankle
1005P	Gel Heel Sleeve L/XL	6 3/4	12	measured at Ankle
1007P	Gel Lace Pad Sleeve One Size	6 3/4	10	measured at Ankle
1002P	Heel and Lace Pad Sleeve S/M	6 3/4	10	measured at Ankle
1003P	Heel and Lace Pad Sleeve L/XL	6 3/4	12	measured at Ankle
317P	Knee/Elbow Sleeve S/M	6 3/4	13	
318P	Knee/Elbow Sleeve L/XL	6 3/4	13	
329P	Knee/Elbow Sleeve XXL	9 1/2	19	
328P	Bamboo Gel Knee Pad Sleeve S/M Length 11 1/2"	9 1/2"	19"	Thigh/ measured 6" above
		8"	17"	Calf/measured 6" below k
325P	Bamboo Gel Knee Pad Sleeve L/XL Length 11 1/2"	11"	22"	Thigh/measured 6" above
		9"	18"	Calf/measured 6" below k
312P	Ultra Thin Dress Bunion 3mmGel Small	6 1/2	12	measured at foot opening
313P	Ultra Thin Dress Bunion 3mmGel M/L	6 1/2	12	measured at foot opening
314P	Ultra Thin Dress Bunion 3mmGel XXL	6 1/2	12 1/4	measured at foot opening
311P	All Gel Bunion - Great Toe	1 1/2"	2 1/2"	
310P	All Gel Bunion - Small Toe	1"	2 3/4"	
303P	Foot Cover - Small	8	10	measured at foot opening
304P	Foot Cover - Medium	8 1/4	10 1/2	measured at foot opening
305P	Foot Cover Large	9	11	measured at foot opening
P/322P/0202	Goalie Shock Pad/Wrist Gel Sleeve S/M	6 3/4	10	measured at palm
P/324P/0204	Goalie Shock Pad/Wrist Gel Sleeve L/XL	6 3/4	12	measured at palm
Product Number	Product Description	Width (in)	Length (in)	Girth/HT (in)
500P	Spreaders Small	1 1/16	7/8	3/8
501P	Spreaders Medium	1 1/8	1	3/8
502P	Spreaders Large	1 3/16	1 1/16	3/8
503P	Separators Small	7/8	1 1/2	1/8
504P	Separators Large	1	2	1/8
328P	Gel Self Adhesive Oval w/Fabric Cover			
205P	Gel Pad Rectangular w/Hook backing	1 3/4	2 1/2	~1/4" thick
206P	Gel Pad Circle w/Hook backing		2 1/2 (dia)	~1/4" thick

85% Nylon; 15% Spandex; Gel Pad: TPE (Thermo Polymer Elastomer)

43% Nano Bamboo Charcoal;17% Lycra;21% Spandex;19% Nylon

Fabric Content: 83% Nylon, 17% Spandex, Gel Pad, TPE (Thermo Polymer Elastomer)

Fabric Content: 96% Nylon; 4% Spandex; Gel Pad: TPE

... Polyester; Gel Pad:TPE (Thermo Polymer Urethane)

973-296-9994

UnlimitedMotionllc@gmail.com
ndudley@unlimitedmotionllc.com